

# Hydropower and Ocean Energy

# 1.2

2E a division of



Introduction		Hydropower		Ocean Energy	
	<b>Subject Areas</b> Hydropower and Ocean Energy	044		<b>Basic Knowledge</b> Hydropower	046
	<b>HM 150.19</b> Operating Principle of a Pelton Turbine	048		<b>HM 150.20</b> Operating Principle of a Francis Turbine	049
	<b>HM 450.01</b> Pelton Turbine	050		<b>HM 450.02</b> Francis Turbine	050
	<b>HM 450C</b> Characteristic Variables of Hydraulic Turbomachines	051		<b>HM 421</b> Kaplan Turbine Trainer	052
	<b>HM 365.31</b> Pelton and Francis Turbine	054		<b>HM 430C</b> Francis Turbine Trainer	056
	<b>ET 270</b> Wave Energy Converter	060		<b>Basic Knowledge</b> Wave Energy	058


 Visit our website at:  
[www.gunt2e.de](http://www.gunt2e.de)




# Subject Areas

# Hydropower and Ocean Energy

Subject Areas

Products

Water's natural flowing movements, such as in rivers and reservoirs, can be used in the production of electricity. Furthermore, both the tidal range (the periodic rise and fall of the sea level) and the energy contained in flow and waves can be used in the ocean sector.

Both types of energy conversion are classed as renewable energies. While the typical use of hydropower has been widespread for hundreds of years, using the ocean for energy is in its infancy.

As the table opposite shows, different learning objectives from turbine engineering can be differentiated in the expanded field of hydropower and ocean energy. The corresponding 2E product is listed in the next column.



Additional trainers in the fields of turbines and fluid mechanics in particular can be found in GUNT catalogue 4 "Fluid Mechanics".

## Hydropower

**Fundamentals of fluid mechanics:**  
energy conversion in water turbines

**HM 150.19**  
Operating Principle of a Pelton Turbine

**HM 150.20**  
Operating Principle of a Francis Turbine

**Turbomachines:**  
measurements on turbines and pumps

**HM 450.01**  
Pelton Turbine

**HM 450.02**  
Francis Turbine

**HM 450C**  
Characteristic Variables of Hydraulic Turbomachines

**Turbines in run-of-river power stations:**  
river and tidal power stations

**HM 421**  
Kaplan Turbine Trainer

Comparing turbine types – asynchronous generator drive

**HM 365.31**  
Pelton and Francis Turbine

How guide vanes affect characteristic turbine curves

**HM 430C**  
Francis Turbine Trainer

## Ocean Energy

How the guide vanes affect characteristic turbine curves

**ET 270**  
Wave Energy Converter

# Basic Knowledge Hydropower



Traditional hydropower systems have been in use for hundreds of years as a source of energy for a wide variety of mechanical applications. As such, hydropower represents a renewable energy source that has been successfully used for a long time. Since the beginning of hydropower-generated electricity, the percentage of electrical energy generated in this way has grown to around one quarter of all the electricity used worldwide.

However as the number of turbines in use increases, and with it the necessary retaining dams, the significant drawbacks in the overall ecological balance of this technology has become apparent to some extent. Due to geological conditions, some countries such as Norway (99%), Congo (97%) and Brazil (96%)

are able to cover very large proportions of their electrical energy demand with hydropower. By comparison, in Germany only 4% is covered. Brazil is currently home to the world's largest hydroelectric power station, where 18 turbines generate a total output of 12,600 megawatts.




## Decentralised power supply through small hydroelectric power stations


In regions without a central power supply, decentralised, small hydroelectric power stations with an output up to about 5kW offer the possibility of supporting sustainable development in an appropriate manner.

In addition to the typical characteristic variables such as head and flow rate, other aspects such as maintenance issues and accessibility of the installation site are also important to consider when selecting the type of turbine. At heads of 150m and more Pelton turbines are mostly used. At lower heads on the other hand, Francis or Kaplan turbines are preferred.

## Expansion potential for use of hydropower




## Turbine types in hydroelectric power stations



### Pelton turbine

In the Pelton turbine the water "shoots" out of one or more nozzles onto the vanes of the impeller.


Head: 150-2000m  
Flow rate: 0,02-70m<sup>3</sup>/s  
Storage power stations



### Francis turbine

The Francis turbine operates with positive pressure. The guide vanes can be adjusted.

Head: 20-700m  
Flow rate: 0,3-1000m<sup>3</sup>/s  
Dams



### Kaplan turbine


The Kaplan turbine also operates with positive pressure. In this case, guide vanes and impeller vanes can be adjusted.

Head: 2-60m  
Flow rate: 4-2000m<sup>3</sup>/s  
Rivers


1 impeller 2 guide vanes 3 water inlet 4 water outlet

# HM 150.19 Operating Principle of a Pelton Turbine

In order to successfully teach the fundamentals of energy conversion in hydroelectric power stations, an understanding of the fluid mechanical processes directly upstream of the turbine is especially important.



Observe the water flow out of the nozzle onto the impeller through the transparent front panel on HM 150.19. By adjusting the nozzle needle, you can vary the nozzle cross-section and study its effect on the operating behaviour.



**Learning objectives**

- design and function of a Pelton turbine
- determination of torque, power and efficiency
- graphical representation of characteristic curves for torque, power and efficiency

HM 150 Base module




Product No. 070.15019  
More details and technical data:  
[gunt.de/static/s4562\\_1.php](http://gunt.de/static/s4562_1.php)




# HM 150.20 Operating Principle of a Francis Turbine

The Francis turbine is a positive pressure turbine, in which conversion of the water's pressure energy into kinetic energy takes place in the guide vanes and in the impeller. The water is directed to the guide vanes via a spiral shaped housing.



You can adjust the alignment of the Francis turbine's guide vanes to optimise performance.

HM 150.19 and HM 150.20 are part of the HM 150 series on fluid mechanics. We recommend using the HM 150 base module for supplying the turbines with water and for measuring the flow rates.



HM 150 Base Module for Experiments in Fluid Mechanics

**Learning objectives**

- design and function of a Francis turbine
- determination of torque, power and efficiency
- graphical representation of characteristic curves for torque, power and efficiency

Product No. 070.15020  
More details and technical data:  
[gunt.de/static/s4572\\_1.php](http://gunt.de/static/s4572_1.php)





# HM 450.01 HM 450.02

## Pelton Turbine Francis Turbine

With HM 450.01 and HM 450.02 you can systematically investigate energy conversion in various turbine types and compare typical properties. Used in conjunction with the HM 450C trainer, you have a fully equipped system for education in the field of turbomachines.

The water exits the Francis impeller axially. The intake capacity, and therefore the performance of the turbine, can be adjusted by rotating the guide vanes.



# HM 450C

## Characteristic Variables of Hydraulic Turbomachines




HM 450C contains a closed water circuit to supply the turbines. Pressure, flow rate, torque and rotational speed are detected electronically and displayed digitally.

Learning objectives
<ul style="list-style-type: none"> <li>comparison of turbine types</li> <li>determining the mechanical power</li> <li>influence of nozzle cross-section and/or guide vanes</li> <li>plotting characteristic curves</li> <li>calculating the efficiency</li> </ul>

### The HM 450C software

Measured values can be transferred to a PC via USB. The GUNT software provides you with an extensive range of options for displaying and analysing your data.



Product No.  
070.45001  
More details and technical data:  
[gunt.de/static/s4202\\_1.php](http://gunt.de/static/s4202_1.php)



Product No.  
070.45002  
More details and technical data:  
[gunt.de/static/s4201\\_1.php](http://gunt.de/static/s4201_1.php)





Product No.  
070.4500C  
More details and technical data:  
[gunt.de/static/s4242\\_1.php](http://gunt.de/static/s4242_1.php)




# HM 421 Kaplan Turbine Trainer

Kaplan turbines are used at low hydraulic heads to generate electricity. Low heads occur for example, in run-of-river power stations and tidal power stations.

HM 421 provides you with a trainer equipped with a Kaplan turbine in a closed water circuit and a submersible pump.



Impeller                      Guide vanes



### Software

The software for HM 421 allows the most important variables to be captured:

- flow rate
- head
- speed
- torque
- pressure at the inlet and outlet
- temperature

The following variables can be calculated from the measured values:

- hydraulic power
- mechanical power
- efficiency
- head

### Learning objectives

- measurement of turbine characteristics
- determining power characteristics at different speeds (hydraulic power, mechanical power)
- determining the head
- calculating the turbine efficiency
- how blade position affects power and efficiency

Product No.  
070.42100  
More details and technical data:  
[gunt.de/static/s4917\\_1.php](http://gunt.de/static/s4917_1.php)


# HM 365.31 Pelton and Francis Turbine

The HM 365 modular system allows you to investigate the characteristic operating behaviour of various turbine types. Water is supplied by the specially designed HM 365.32 supply unit. The energy generated by the turbine is transferred to the asynchronous generator in HM 365. Further information about the additional possibilities can be found in the data sheets for the respective devices.



The **GUNT software** provides an intuitive visualisation of current measurement data in a system diagram and enables continuous data acquisition via a USB connection. Of course, graphical representations and calculations are also provided for analysis of the measurement data.

Product No.  
070.36531  
More details and technical data:  
[gunt.de/static/s4233\\_1.php](http://gunt.de/static/s4233_1.php)



HM 365.31: Pelton Turbine



HM 365.31: Francis Turbine



HM 365 Universal Drive and Brake Unit

HM 365.32 Turbine Supply Unit

## Learning objectives

- comparison of impulse and reaction turbines
- determination of mechanical and hydraulic power
- determination of efficiency
- recording of characteristic curves
- influence of Pelton turbine nozzle cross-section on characteristic values
- influence of Francis turbine guide vane position on characteristic values

HM 365 has a three-phase asynchronous motor which is used as a generator. Thanks to the possibility of controlling the load via the speed or torque, the turbine being studied can be operated at the optimum operating point with varying hydraulic power.

HM 365.32 is equipped with sensors for pressure and flow rate. A power pump – inside an enclosed water circuit – simulates the gradient of a hydroelectric power station.

# HM 430C Francis Turbine Trainer

Francis turbines are used for hydropower at medium heads and medium flow rates. HM 430C provides you with the facility for using a DC generator for energy conversion.




Torque measurement

DC generator as load


Francis turbine

Centrifugal pump




### Software


The GUNT software displays torque, speed and pressure to monitor the current system state.



You can easily study the behaviour of the turbine with different guide vane positions.



A force sensor detects the torque at the generator's flexible mounting.



The mechanical power emitted, the hydraulic power supplied and the efficiency can all be calculated from the measurement data. The chart shows how the efficiency depends on the speed for three different guide vane positions.

### Learning objectives

- investigation of conversion of hydraulic into mechanical energy
- determination of torque and speed of turbine shaft
- determination of mechanical and hydraulic power
- determination of efficiency
- recording of characteristic curves
- investigate the effect of guide vane position
- calculate velocity triangles

Product No.  
070.430C0  
More details and technical data:  
[gunt.de/static/s4585\\_1.php](http://gunt.de/static/s4585_1.php)



# Basic Knowledge Wave Energy




The ocean's waves contain an inexhaustible reserve of energy. They are caused by wind, gravitational forces and atmospheric pressure differences.

The International Energy Agency estimates the possible global contribution of wave energy to the power supply at more than 10%.

The main challenge in constructing wave power stations is not least that of designing systems that can withstand the sometimes destructive natural conditions. The integration of chamber systems following the principle of the oscillating water column (OWC) has proven promising in existing coastal defence structures.

Annual average power of ocean waves (kW/m)



- < 5
- 5 – 10
- 10 – 15
- 15 – 20
- 20 – 30
- 30 – 40
- 40 – 60
- > 60

The map shows the average annual wave power. It assumes power along a coastline or along a wave crest. The power densities are specified in kW/m. It should be noted that high powers can be found at latitudes far from the equator and on the west coasts of the continents.


Source: Centre for Renewable and Sustainable Energy Studies, Stellenbosch University

Linear wave theory provides an estimate for the energy flow of a wave:

$$P \sim T * H^2$$


We can see that the power **P** is linearly dependent on the period length **T** and quadratically dependent on the height of the wave **H**.

Distribution of the wave energy supply



When designing systems to use wave energy, the distribution of the wave energy supply is particularly important. Results from global studies show that the largest percentage of wave energy can be assigned to a period length between 1 and 30 seconds.

Basic mechanical principles for using wave power



Systems for using wave energy which have been proposed in the past and also partially put into practice can be classified in the following categories according to the underlying principle:



- 1 floating systems
- 2 chamber systems (OWC)
- 3 overflow systems

# ET 270 Wave Energy Converter

Wave power stations use the energy of ocean waves to generate electricity. Systems designed according to the principle of the oscillating water column (OWC) are fitted with a chamber which opens out below sea level. At the other end is a turbine. The continuous wave motion produces an oscillating water column within the chamber, which causes the air mass above it to move. The airflow generated in this way drives the turbine, which generates electricity via a generator.

In addition to the usual components of an OWC wave power station, ET 270 contains a water channel with a controllable wave generator for laboratory experiments. The GUNT software enables the acquisition and representation of relevant measured values such as wave height and speed of the generator.


Modelled on large-scale power stations, ET 270 also uses a Wells turbine. This type of turbine functions regardless of the flow direction of the air. Consequently, it is possible to use the energy of the airflow as it moves both up and down.



The speed of the Wells turbine can be specified by controlling the generator. As a result, it is possible to adjust the most efficient operating point for energy recovery.

**Learning objectives**

- familiarisation with wave power plants
- understanding of energy generation from wave motion
- measurement of wave motions
- familiarisation with the Wells turbine
- optimisation of operating behaviour



**Software**  
Measurement of the wave height

Product No.  
061.27000  
More details and technical data:  
[gunt.de/static/s5372\\_1.php](http://gunt.de/static/s5372_1.php)



## Energy



## Hydropower and ocean energy

**HM 150.19**  
 Operating principle  
 of a Pelton turbine

Model of an impulse turbine with adjustable nozzle; determination of efficiency

Recommended for water supply: HM 150 Base module for experiments in fluid mechanics

Order No.: 070.15019


**HM 150.20**  
 Operating principle  
 of a Francis turbine

Model of a reaction turbine with adjustable guide vanes and determination of the efficiency

Recommended for water supply: HM 150 Base module for experiments in fluid mechanics

Order No.: 070.15020


**HM 450.01**  
 Pelton turbine


Model of an impulse turbine with speed and torque measurement

Order No.: 070.45001


**HM 450.02**  
 Francis turbine

Model of a reaction turbine with speed and torque measurement; adjustable guide vanes

Order No.: 070.45002


**HM 450C**  
 Characteristic variables of hydraulic turbomachines

Determination of output and efficiency of turbines and pumps; demonstration of a pumped storage plant

Order No.: 070.450C0


**HM 421**  
 Kaplan turbine  
 trainer

Four-bladed Kaplan turbine with guide vane adjustment for varying power

Order No.: 070.42100


**HM 430C**  
 Francis turbine  
 trainer

Characteristics of a powerful Francis turbine with adjustable guide vanes

Order No.: 070.430C0


**HM 365.31**  
 Pelton and Francis  
 turbine

Comparison of impulse and reaction turbines

Order No.: 070.36531


**HM 365.32**  
 Turbine supply unit

Water supply for HM 365.31

Order No.: 070.36532


**HM 365**  
 Universal drive and  
 brake unit

Core component for experiments on various driving and driven machines

Order No.: 070.36500



Trainer for turbines with Pelton turbine HM 365.31, supply unit HM 365.32 and brake unit HM 365

**ET 270**  
 Wave energy converter

Turbine unit with Wells turbine and electric generator; configurable wave generator


Order No.: 061.27000



# Hydropower

## HM 450C Characteristic variables of hydraulic turbomachines

The trainer HM 450C is able to operate pump and turbine at the same time. Relevant measured values are recorded contemporaneously at both turbomachines. Thus the trainer can be used as a **pumped storage plant**.



HM 450.01  
Pelton turbine

HM 450.02  
Francis turbine

HM 450C  
Trainer with  
water turbines



Needle nozzle and  
impeller of the Pelton  
turbine



Adjusting knob for the needle  
nozzle



Position of the guide vanes in the  
Francis turbine



Vanes and impeller of the  
Francis turbine

## HM 288 Experiments with a reaction turbine



## HM 289 Experiments with a Pelton turbine



## HM 291 Experiments with an action turbine



HM 290  
Base unit for  
turbines

## HM 287 Experiments with an axial turbine

